



MIPAR™

Release Notes

v1.0.2

Bugs Fixed

Image Processor

- Fixed bug with optimizing *Retain Dilation* Recipe steps
- Fixed bug with “playing” *Area Fraction* Recipe steps
- Fixed bug with display of several Recipe steps such as *Watershed* and *Clean Boundaries*
- Fixed bug with loading Recipes which had an *Intercepts* measurement step whose “Save All Intercepts” parameter was set to “Yes”

3D Toolbox

- Fixed bug with *Separate Features 3D*
- Fixed bug when using B/W reference images
- Fixed bug with *Remove Specific Features*

Other Changes

- Default mode for 3D zoom is now “Camera” for improved behavior
- Import menu is now hidden in the 3D Toolbox