




MIPAR™

Release Notes

v1.3.1

Auto Mag Calibration



Setting your pixel size has never been easier. MIPAR now automatically detects the scale bar for you!


Auto → **Position the box** → **Confirm**

Type the known distance and you're done!

Find a scale bar it doesn't detect? Let us know and we'll keep refining!

[Watch Tutorial »](#)

Integrated License Activation



MIPAR licenses of all types can now be activated directly within MIPAR from the new Activation Window!

Appears after launch if no active license.


You may also access it in any app from Help > Activate License

Bulk Measurement Stats



Processed
Image Editor

FEATURE
MEASURES



MEANS per Image	

STATS per Image	

- ✓ Mean
- ✓ Min
- ✓ Max
- ✓ StdDev

Now when generating feature measurements from hundreds of images in the Processed Image Editor, in addition to the spreadsheet of feature measurements for all images, you'll see a "stats" sheet for all images.

Fixes and Changes

- Minor bug fixes throughout
- "Adjustment Connectivity" setting removed for simplicity in *Watershed* and *Separate Features*
- *Crop Image* steps early in the Recipe now effect subsequent *Channel Operation*, *Restore Original*, and *Color Cluster* steps
- Moved *Replace With* and *Mark Center* functions to "Clean-Up" menu